

Curera[®] Back Pillow

Pillow with wing for lateral positioning/side support



 **curera**[™]
positionering

 care of sweden
SUPPORTING LIFE

Curera® Back Pillow

If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.


Curera® Positioning is a range of stable pillows, rolls, and wedges that meet most needs for support and relief. They are also a great help when changing dressings and for patient care.

Curera® Back Pillow is a winged-pillow for back relief. The pillow is placed behind the back, and the wing rests under the user to set the position.

Working with Curera® Positioning is simple. Contact Care of Sweden for training and further information.



POSITIONING EXAMPLE

CLEANING	FACTS
<p>Inner pillow: Machine wash at max. 95 °C. Can be tumble-dried or dried in an airing cupboard; avoid over-drying. For faster drying, run the spin cycle three or four times.</p>	<ul style="list-style-type: none">• Filling: 50 % HR foam/ 50 % polyester fibre• Warranty: 1 year• Size: 50 x 85 cm <p>CE</p> <p> Always read the instruction for use prior to usage of the product.</p>

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2019. All rights reserved. Reproduction prohibited. Last updated: 2021-06-01