


Curera® Torso Support Pillow

Pillow for torso support in seated position



 **curera**[™]
positionering

 **care of sweden**
SUPPORTING LIFE

Curera® Torso Support Pillow

If a person makes too few spontaneous movements, there is an increased risk of pressure ulcers. Regular changes of position reduce the risk of pressure ulcers. Position as required in both reclining and seated position.

Curera® Positioning is a range of stable pillows, rolls, and wedges that meet most needs for support and relief. They are also a great help when changing dressings and for patient care.

Our hygiene covers and cotton covers are also available as accessories to the pillows and rolls. The hygiene cover does not allow liquid or air to penetrate, but is vapour permeable (vapour is wicked away).


Make sure that the patient is positioned correctly so that airways are unobstructed.

Curera® Torso Support Pillow is used to support the torso in the seated position and provides good underarm support. The pillow can be placed from the front and from behind.

Working with Curera® Positioning is simple. Contact Care of Sweden for training and further information.



POSITIONING EXAMPLE

CLEANING	FACTS
<p>The hygiene cover can be wiped clean with detergent and/or disinfectant. Machine wash at max. 95 °C.</p> <p>Inner pillow: Machine wash at max. 95 °C. Can be tumble-dried or dried in an airing cupboard; avoid over-drying. For faster drying, run the spin cycle three or four times.</p> <p>Cotton cover: Machine wash at max. 70 °C.</p>	<ul style="list-style-type: none"> • Sizes: 85x65x30 cm • Filling: 50 % HR foam/50 % polyester fibre • Cover: Polyurethane/polyester. • Accessories: Cotton cover. • Warranty: 1 year • Accessories: Hygiene cover: Polyurethane/polyester Cotton cover <p style="text-align: right;">  Always read the instruction for use prior to usage of the product. </p>

Care of Sweden and Curera are registered trademarks of Care of Sweden AB. © Care of Sweden AB, 2013. All rights reserved. Reproduction prohibited. Last updated: 09/04/2018